

BARISELECT®

Low Air Loss Mattress Replacement



User-Friendly Control Unit

Rapid inflation, easy set-up and a motor so quiet, the patient will hardly know it's there!

The SelectAir® System Quietly circulates air through individual cells and allows air to escape from the mattress.

Wicking Action: Air escapes through tiny holes on the top of the mattress and circulates under a vapor-permeable cover to wick moisture away from patient.

BARISELECT® Features

Two extra inches of mattress depth (8"), extended widths and specialized software accommodate extra weight capacity.

Pulsation (alternating pressure): Provides two cycle settings that continually readjust the system to ensure patient comfort.

Manual Mode allows you to input the patient's actual weight for optimal immersion.

LCD Screen with back-lit keys allowing system features to be visible in dark or dimly lit environments.

Auto Comfort: Takes the guess work out of programming. During set-up, mattress sets optimal pressure according to patient's weight distribution.

Zoned, "true" low air loss used in acute care hospitals nationwide

Weight Limit: 1,000 lbs. (453 kg)

Mattress: Fits most bariatric bed frames. Custom sizes available. Inquire for details.

Cover: Water resistant, vapor-permeable, flame resistant and biocompatible

Warranty: One Year Limited Warranty on control unit and soft goods

Medicare Code: E0277 Group 2 Support Surface

LIT-0005 | R100121

Manufactured in the USA by
MOXI™
ENTERPRISES
Innovative health solutions

1939 S. Vandeventer Ave.
St. Louis, MO 63110

T 1-800-942-7678
E info@moxiusa.com
W moxiusa.com